95x210mm

Tolaccea

Hot & Cold Therapy Compression Ice Foot Sleeve

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WARNING:
To avoid danger of suffocation, keep this away from children and babies!

[EC] [REP] Company Name: REPLAB SARL
Address and Post Code: 99 AVENUE ACHILLE
PERETTI 92200 NEUILLY-SUR-SEINE France
Contact: XIAOLAN LIU
Phone: +33767189223
Email: EUREP@REPLAB.CC

[UK] [REP] Company Name: REPLAB CONSULTING LTD
Address and Post Code: 376 NEWPORT ROAD CARDIFF,
WALES, CF23 9AE
Contact: HUIHUA LU
Phone: +447395577022
Email: UKREP@REPLAB.CC

UK/US/CA

FUNCTIONS

1. Tolaccea hot & Cold therapy ice foot sleeve provides 360° coverage;
2. The flexible ice foot sleeve is in ergonomic designs perfect for foot and it is wearable without straps;
3. This ice foot sleeve is reusable, also microwave and freezer safe, allowing you to use hot or cold therapy depending on your needs;
4. The ice foot sleeve is made of non-toxic material, latex-free, and leak proof, with no smells after use.
5. It can help to relief pain for : achilles tendonitis, plantar fasciitis, sprained ankles or feet

OPERATION INSTRUCTION

- 1.Cold therapy**
- Place the foot sleeve flat in a refrigerator at -20 °C for 1-2 hours for immediate use. If colder temperature needed, frozen it longer with lower temperature.
 - Take out the foot sleeve from freezer, wipe off the moisture and pull it over the required area.
 - Wear foot sleeve and keep in place for approximately 15-20 minutes
- 2.Hot therapy**
- Start with room temperature pack, place the foot sleeve flat in microwave with a clean plate.
 - Microwave the foot sleeve for 10 seconds, if more heat is needed microwave in 5 second increments.
 - Wear it and keep in place for approximately 15-20 minutes

HOW TO CLEAN IT

Hand wash with damp towel. Do not submerge sleeve in water. Lay flat and air dry, do not machine wash.

WARNING!

- 1.Limit cold/hot therapy to 30 minutes to prevent frostbite or burns. If not cold enough, refreeze at -20°C.
- 2.Do not microwave over 45 seconds. Check temperature before use; cool if too hot to avoid burns.
- 3.Repeated heating may dry out the gel, reducing lifespan. Monitor while microwaving to avoid overheating. Air dry if damp—never use a dryer.
- 4.For external use only. Avoid areas with poor circulation, open wounds, or infections unless advised by a professional.
- 5.If gel contacts skin, wash thoroughly. Seek medical advice for injuries, and discontinue if discomfort persists.
- 6.Store in therapy wrap when not in use to prevent condensation. Ensure dry before storage.
- 7.Keep away from children and pets.

AFTER SALE

Tolaccea guarantees 1 year warranty for the quality problem. Any questions please freely email us: service@tolaccea.com

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