21*9.5cm

Tolaccea

Tightened Full Leg

Ice Pack Comes with Expert Cooling Gel and Velcro

PRODUCT MANUAL

HOT PRODUCTS

Elbow & Knee Compression Sleeve Ice Pack
Colour: Black, Blue
Size: S, M, L, XL; Unisex

Ankle Foot Ice Pack Wrap
Colour: Black, Blue
Size: L, XL; Unisex

Stress Relief Balls
Color: Black, Blue, Light Blue
Size: 2.4in

3D SHOULDER ICE PACK COOLING GEL BRACE
Colour: Black
Size: Universal , L Unisex

Half Finger Arthritis Compression Ice Glove
Colour: Black
Size: S , M; Unisex

Reusable Gel Cold Therapy Sock
Colour: Black
Size: M, L; Unisex

Thumb Wrist Ice Pack with Adjustable Wrist Strap
Cour: Black
Size: One size fits all; Unisex

Foot Ice Packs Big Toe for Pain Relief
Colour: Black
Size: One size fits all; Unisex

Wrist Ice Pack Wrap
Colour: Black
Size: S , M; Unisex

Foot Ice Pack for Toes
Colour: Black
Size: One size fits all; Unisex

- EN -

FUNCTIONS

Full-Leg Coverage: 360° wrap-around design for effective icing of the entire leg muscles.
Multi-Region Adaptability: Flexible use for knees, ankles, calves, thighs, and feet.
Post-Workout Recovery: Designed for relaxation and muscle recovery after exercise, relieving fatigue.
Comfortable Fit: Ergonomic design ensures a snug fit for maximum comfort and effectiveness.
Efficient Cooling: Provides rapid cooling and maintains a consistent low temperature, ideal for muscle recovery.

OPERATIONINSTRUCTION

Cold therapy

- Place the foot sleeve flat in the freezer at -20°C (0°F) for 1-2 hours for immediate use. If you need it colder, freeze it longer or at a lower temperature.
- Remove the foot sleeve from the freezer, wipe off any excess moisture, and slide it onto the desired area.
- Keep the foot sleeve in place for approximately 15-20 seconds.

Hot therapy

- Remove the product from the packaging.
- Fold it according to the marked sections on the packaging. Heat each side for 40 seconds in the microwave.
- If additional heat is required, unfold the product, refold it inside-out, and heat for another 10-20 seconds.
- Ensure the microwave is set to a standard power level of 700W or below.
- Wear the foot sleeve and keep it in place for approximately 15-20 minutes.

HOW TO CLEAN IT

To clean the sleeve, hand wash it with a damp towel and avoid submerging it in water. Lay it flat to air dry, and do not machine wash or dry.

WARNING!

- Limit cold/hot therapy to 30 minutes to prevent frostbite or burns. If not cold enough, refreeze at -20°C.
- Do not microwave over 45 seconds. Check temperature before use; cool if too hot to avoid burns.
- Repeated heating may dry out the gel, reducing lifespan. Monitor while microwaving to avoid overheating. Air dry if damp—never use a dryer.
- For external use only. Avoid areas with poor circulation, open wounds, or infections unless advised by a professional.
- If gel contacts skin, wash thoroughly. Seek medical advice for injuries, and discontinue if discomfort persists.
- Store in therapy wrap when not in use to prevent condensation. Ensure dry before storage.
- Keep away from children and pets.

AFTER SALE

Tolaccea guarantees 1 year warranty for the quality problem.
Any questions please freely email us: service@tolaccea.com.

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Company Name : REPLAB CONSULTING LTD