96mm

EN-1

FUNCTIONS

- 1. Immediate Relief: Easy to use cold therapy sock for foot care. Ideal for sore feet, Arthritis, Achilles Tendonitis, Plantar Fasciitis, Sprained, Swelling, and Post-Surgery Recovery. An excellent at-home cold therapy solution that provides relief in minutes.
- 2. Easy To Use: No need to replace/fill ice packs. Say no to messy and inconvenient methods like a zipped bag of lumpy ice which never seemed to cool the right spot. It remains soft after freezing and does not affect free movement after wearing. Experience our cool socks for long-lasting relief.
- 3. Hot & Cold Therapy: Tolaccea Cold therapy sock for cold and heat therapy depending on your needs. The expert gel locks in cold or heat for a long time for better penetration into your foot and ankle.
- 4. Professional Cutting-Edge Materials: It is made of Lycra fabric and expert-grade thick solid gel, durable, soft, skin-friendly, and reusable. It can effectively relieve foot inflammation, pain, and postoperative recovery.
- 5. Certified Safe and Fully Guaranteed: Using a hot or cold compress is a natural and effective way to fight inflammation and promote faster recovery. Moreover, Tolaccea's customer support team provides after-sales support and services to customers worldwide. If you have any questions, please feel free to contact us, we will give you a satisfactory solution.

OPERATION INSTRUCTION

1.COLD THERAPY

- Place the sock flat in a refrigerator at -20°C for 4-6 hours for immediate use. If colder temperature needed, frozen it longer with lower temperature;
- Take out from freezer, wipe off the moisture and pull it over the required area. Wear and keep it in place for approximately 15-20 minutes.

2.HOT THERAPY

- With the support at room temperature, place the sock in center of microwave.
- Microwave for 50 seconds and test the sock temperature for comfort.
- If the sock is not warm enough, reheat for additional 10 seconds.DO NOT heat for more than 90 seconds. The sock should be warm, not hot. If sock is too hot for comfort, allow it to cool until a comfortable temperature is reached.

CAUTION:

Always monitor the sock while heating in microwave. The sock will absorb water. If it becomes damp, allow to air dry. Do not place the sock in dryer. Do not apply to infected skin or open wounds without medical consultation. If gel contacts body, wash skin thoroughly with soap and water.

EN-2

WARM REMINDER:

According to the doctor's suggestion, the best time to use cold/heat therapy is about 15–20 minutes. Too cold, overheating, and prolonged cold therapy or heat therapy may cause damage to your body, please pay attention to controlling the temperature of use of the product and the duration of use.

CARE INSTRUCTIONS:

Use damp cloth to wipe down. Lay flat and air dry. Do not wash, bleach, dry clean or iron.

⚠ WARNING!

- 1.Limit cold/hot therapy to 30 minutes to prevent frostbite or burns. If not cold enough, refreeze at -20°C.
- 2.Do not microwave over 45 seconds. Check temperature before use; cool if too hot to avoid burns.
- 3.Repeated heating may dry out the gel, reducing lifespan. Monitor while microwaving to avoid overheating. Air dry if damp—never use a dryer.
- 4.For external use only. Avoid areas with poor circulation, open wounds, or infections unless advised by a professional.
- 5.If gel contacts skin, wash thoroughly. Seek medical advice for injuries, and discontinue if discomfort persists.
- 6.Store in therapy wrap when not in use to prevent condensation. Ensure dry before storage.
- 7.Keep away from children and pets.

AFTER SALE

Tolaccea guarantees 1 year warranty for the quality problem. Any questions please freely email us:service@tolaccea.com

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