100mm

100mm

Tolaccea**Toes and Foot Ice Sleeve**

Operation Instruction

@tolaccea

www.tolaccea.com

#tolaccea

service@tolaccea.com

Company Name: REPLAB SARL

Address and Post Code: 99 AVENUE ACHILLE
PERETTI 92200 NEUILLY-SUR-SEINE France

Contact: XIAOLAN LIU

Phone: +33767189223

Email: EUREP@REPLAB.CC

Company Name: REPLAB CONSULTING LTD

Address and Post Code: 376 NEWPORT ROAD CARDIFF,
WALES, CF23 9AE

Contact: HUIHUA LU

Phone: +447395577022

Email: UKREP@REPLAB.CC

quelairedemesdechets.fr

EN

Functions:

- Tolaccea toes and foot ice pack provides 360° coverage;
- The flexible ice foot is in ergonomic designs perfect for feet and it is wearable without straps;
- This ice foot sleeve is reusable, also microwave and freezer safe, allowing you to use hot or cold therapy depending on your needs;
- The ice pack is made of non-toxic material, latex-free, and leak proof, with no smells after use.
- It can help to relief pain for : toes and foot

Operation Instruction:**•Cold therapy**

Place the ice pack flat in a refrigerator at -20°C for 4-6 hours for immediate use. If colder temperature needed, frozen it longer with lower temperature; Take out the ice pack from freezer, wipe off the moisture and pull it over the required area. Wear and keep it in place for approximately 15-20 minutes.

•Hot therapy

Start with room temperature, place the ice pack flat in microwave with a clean plate; Microwave it for 50 seconds, if more heat is needed microwave in 5 second increments; Wear and keep it in place for approximately 15-20 minutes.

How to clean it

Hand wash with damp towel. Do not submerge sleeve in water. Lay flat and air dry, do not machine wash.

-1-

WARNING!

- 1.Limit cold/hot therapy to 30 minutes to prevent frostbite or burns. If not cold enough, refreeze at -20°C.
- 2.Do not microwave over 45 seconds. Check temperature before use; cool if too hot to avoid burns.
- 3.Repeated heating may dry out the gel, reducing lifespan. Monitor while microwaving to avoid overheating. Air dry if damp—never use a dryer.
- 4.For external use only. Avoid areas with poor circulation, open wounds, or infections unless advised by a professional.
- 5.If gel contacts skin, wash thoroughly. Seek medical advice for injuries, and discontinue if discomfort persists.
- 6.Store in therapy wrap when not in use to prevent condensation. Ensure dry before storage.
- 7.Keep away from children and pets.

After Sale

Tolaccea guarantees 1 year warranty for the quality problem. Any questions please freely email us: service@tolaccea.com

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