21\*9.5cm

Tolaccea

3D DOUBLE COLDFOCUS

GEL COLD EYE MASK

PRODUCT MANUAL

## HOT PRODUCTS

**Elbow & Knee Compression Sleeve Ice Pack**  
Colour: Black, Blue  
Size: S, M, L, XL; Unisex

**Ankle Foot Ice Pack Wrap**  
Colour: Black, Blue  
Size: L, XL; Unisex

**Full Leg Ice Pack**  
Colour: Black;  
SIZE: Universal;L;  
Unisex

**3D SHOULDER ICE PACK COOLING GEL BRACE**  
Colour: Black  
Size: Universal, L  
Unisex

**Half Finger Arthritis Compression Ice Glove**  
Colour: Black  
Size:S, M; Unisex

**Reusable Gel Cold Therapy Sock**  
Colour: Black  
Size: M, L; Unisex

**Thumb Wrist Ice Pack with Adjustable Wrist Strap**  
Colour: Black  
Size: One size fits all; Unisex

**Foot Ice Packs Big Toe for Pain Relief**  
Colour: Black  
Size: One size fits all; Unisex

**Wrist Ice Pack Wrap**  
Colour: Black  
Size: S, M; Unisex

**Foot Ice Pack for Toes**  
Colour: Black  
Size: One size fits all; Unisex

**3D Ankle Ice Pack**  
Colour: Black PRO;  
SIZE: L,XL;  
Unisex

**3D Elbow&Knee Ice Pack**  
Colour: Black PRO;  
SIZE: S,M,L,XL;  
Unisex

- EN -

### DESCRIPTION

Experience soothing eye relief with 3D enhanced cooling gel design, delivering deep, long-lasting cold therapy. For intensified cooling and comfort, frozen gel packs to increase compression and maximize the chill. Ideal for relieving eye fatigue, puffiness, dark circles, migraines, and post-surgery recovery.

### HOT THERAPY INSTRUCTIONS

- Microwave the product for 20-25 seconds at 700 watts. Test for comfort.
- If additional heat is needed, flip the product and heat in 10-second intervals, up to 1 minute per side.

### COLD THERAPY INSTRUCTIONS

- Freeze the 3D eye mask at -20°C for 1-2 hours.
- For enhanced cooling, insert the included liquid gel packs.
- Wipe off any moisture before applying for 10 minutes or as needed.

### CARE INSTRUCTION

Use damp cloth to wipe down. Lay flat and air dry. Do not wash, bleach, dry clean or iron.

### WHAT IN BOX

3D Eye Mask\*1  
Eye Mask Product Instruction\*1

### WARNING!

1. Limit cold/hot therapy to 30 minutes to prevent frostbite or burns. If not cold enough, refreeze at -20°C.
2. Do not microwave over 45 seconds. Check temperature before use; cool if too hot to avoid burns.
3. Repeated heating may dry out the gel, reducing lifespan. Monitor while microwaving to avoid overheating. Air dry if damp—never use a dryer.
4. For external use only. Avoid areas with poor circulation, open wounds, or infections unless advised by a professional.
5. If gel contacts skin, wash thoroughly. Seek medical advice for injuries, and discontinue if discomfort persists.
6. Store in therapy wrap when not in use to prevent condensation. Ensure dry before storage.
7. Keep away from children and pets.

### AFTER SALE

Tolaccea guarantees 1 year warranty for the quality problem. Any questions please freely email us: service@tolaccea.com.

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Company Name: REPLAB CONSULTING LTD